



Endo tip for broken instrument removal,DTE ED15,Satelec ET20,NSK E7 type

Short endodontic tip for root canal retreatment, indicated for removing broken instruments inside the root canal. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG128

Categories: Satelec - DTE - NSK

Description

Short endo tip for scaler, for root canal retreatment, DTE ED15, Satelec ET20 and NSK E7 compatible. It is indicated for removing fractured instruments inside the root canal.

Specifications

- Compatible with DTE ED15, Satelec ET20 and NSK E7
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use