



Long endo retreatment tip compatible with DTE ED4, Satelec ET40

Long endodontic tip for root canal retreatment, indicated for removing obstructions and root filling material from the root canal.

Available

SKU: PG118

Categories: Satelec - DTE - NSK

Description

Long endo tip for scaler, for root canal retreatment, DTE ED4 and Satelec ET40 compatible. It is indicated for removing obstructions and root filling material residues from the root canal while supporting controlled access during retreatment.

Specifications

- Compatible with DTE ED4 and Satelec ET40
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use