



Root canal filling removal endo tip compatible with DTE ED5, NSK E5

Endo tip for root canal retreatment, indicated for removing obstructions and root filling material from the root canal. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG119

Categories: Satelec - DTE - NSK

Description

Endo tip for scaler, for root canal retreatment, DTE ED5 and NSK E5 compatible. It is indicated for removing root canal fillings, obstructions and filling material residues from the root canal.

Specifications

- Compatible with DTE ED5 and NSK E5
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use