


Endo tip for foreign object removal, DTE ED6, Satelec ET25, NSK E6 type



Endo tip for scaler, used to access areas not reached by other tips, Long Canals, SATELEC, DTE, NSK compatible, autoclavable 134°C.

Available

SKU: PG121

Categories: Satelec - DTE - NSK

Description

Thin, flexible endo tip for scaler, for root canal retreatment, DTE ED6, Satelec ET25 and NSK E6 compatible. It is indicated for removing foreign objects from deep root canals while supporting controlled access during treatment.

Specifications

- Compatible with DTE ED6, Satelec ET25 and NSK E6
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use