



## Short endo retreatment tip compatible with DTE ED7

Short, flexible endodontic tip for root canal retreatment, indicated for removing foreign objects from the root canal.

**Available**

**SKU:** PG120

**Categories:** Satelec - DTE - NSK

### Description

---

Short, flexible endo tip for scaler, for root canal retreatment, DTE ED7 compatible. It is indicated for removing foreign objects from the root canal, including areas where controlled access is required.

### Specifications

- Compatible with DTE ED7
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use