



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/instrument-disinfectant-metasys-id-n-refill/>



Instrument disinfectant Metasys ID N, refill

Instrument disinfectant Metasys ID N, refill, 1 x 5 l canister

133,94 €

Available on request

SKU: 60030034

Categories: Instruments

Description

Instrument disinfectant Metasys ID N refill.

A concentrate for the preventive cleaning and disinfection of invasive and non-invasive dental instruments. One 5-liter bottle, diluted at 1%, yields 500 liters of ready-to-use solution. Free from aldehydes and phenols, this product is listed in the VAH directories and offers a diluted solution lifespan of 14 days.

Composition and Specifications:

- **Advanced Active Ingredients:** Quaternary ammonium compounds, alkylamine, and corrosion inhibitors.
- **Certifications:** CE 0044 compliant.

Packaging Instrument disinfectant Metasys ID N refill: 1 x 5 l canister.