



Head Tracking Mouse with voice dictation, wireless, single pedal

Intera System is a wearable mouse for dentists. Control your Windows computer with a motion-tracking module, a single pedal and voice dictation. The cursor follows the movement of the module, reducing reliance on the mouse and keyboard during clinical work. Includes two microphones. Wireless range: up to 10 m.

Shipping expected in approximately 10 days

SKU: FK802MTR

Categories: All-in-One Medical PCs

Description

Intera System for Windows with a single pedal is Tecnomed Italia's wireless **wearable mouse for dentists**, for controlling the computer during clinical work in the dental practice. Move the cursor with the wearable tracking module, use the foot controls to interact with the computer and dictate clinical notes instead of typing.

Position the tracking module to suit your movement

Position the tracking module to suit the movement you find most comfortable. The module translates your movement into cursor control and can be worn in different positions on the body.

A custom ring mount worn under the glove

Dr Barbara Sabiu used her own 3D printer to create a ring mount that holds the module on her finger under the glove. In this configuration, hand movement controls the cursor. This is one example of adapting the module's position to suit the user's preferred movement. The ring mount is Dr Sabiu's own creation and is not supplied with Intera System.

Wearable cursor control and foot controls

The FK802MTR version is compatible with Windows. The single pedal works alongside the wearable pointer and activates speech recognition. The wearable module moves the cursor in response to the user's chosen movement.

Use the tracking module and foot control to view case images, select documents and browse folders. Dictate notes using the voice input function.



Voice dictation for clinical notes

Two microphones are included for voice input. Activate speech recognition with the foot pedal to dictate information into the patient's clinical record while retaining cursor control through the tracking module.

Less reliance on the mouse and keyboard

Intera System combines motion tracking, foot controls and voice input to access digital content and enter notes, reducing the need to switch to the mouse and keyboard during computer use.

Windows compatibility and wireless connection

This configuration is compatible with Windows and offers a wireless range of up to 10 m. Product code FK802MTR identifies the version with a single pedal. Choose the version that matches your workstation's operating system and the required foot control configuration.

Intera System technical specifications

Product code	FK802MTR
Type	Computer control system with motion tracking, foot controls and voice dictation
Compatibility	Windows PC
Connection	Wireless
Cursor control	Movement of the wearable module, positioned to suit the chosen movement
Foot control	Single pedal
Voice dictation	Included, activated by foot pedal
Microphones included	2
Wireless range	Up to 10 m