



Head Tracking Mouse with voice dictation, wireless, PC Win compatible

Wireless mouse that allows you to perform all the operations of a traditional mouse using head movements and a foot pedal. In addition, the voice dictation function.

Shipping expected in approximately 10 days

SKU: FK801MTR

Categories: All-in-One Medical PCs

Description

Head Tracking Mouse Wireless Windows compatible: Revolutionizing Surgical Interaction

Tecnomed Italia introduces the Head Tracking Mouse Wireless, an innovative solution that transforms PC interaction during surgical procedures, ensuring hygiene and ease of use without compromise.

Utilizing advanced technology, this device tracks head movements to control the cursor hands-free. Foot pedals replicate mouse clicks, providing a safe and intuitive experience.

Key Features and Strengths Head Tracking Mouse:

- **Hygiene and Safety:** By using head tracking and foot pedals, the device minimizes contact with potentially contaminated surfaces, maintaining a sterile and safe operating environment.
- **Operational Freedom:** Navigate files and folders on the PC without using your hands, keeping focus on the procedure and enhancing the concentration of the medical professional.
- **Voice Dictation:** The two included microphones enable voice recognition activated by foot pedals, streamlining clinical documentation directly into the patient's medical record without interruptions.
- **Improved Efficiency:** The system reduces clinical annotation time by up to 80%, optimizing workflow and real-time management of clinical data.
- **Versatility:** Facilitates access to imaging files, captures photos during procedures, and simplifies documentation for assistants, enhancing the entire operational process.
- **Document Precision:** Enhances the quality of clinical notes, contributing to accurate patient record management and ensuring information is always up-to-date.

Technical Specifications Head Tracking Mouse:

- **Connection:** Wireless with head movement recognition, providing cable-free freedom of movement.



- **Pedals:** Three pedals for right click, left click, and voice recognition activation, simplifying digital interaction.
- **Compatibility:** Compatible with Windows PCs, adaptable to various configurations.
- **Equipment:** Includes two microphones for voice commands and dictation during use, enhancing convenience.
- **Range:** Wireless connection up to 10 meters, ideal for operating rooms of various sizes.

The **Head Tracking Mouse Wireless Windows Compatible** offers an advanced solution that combines hygiene and innovation, enhancing digital interaction in modern, dynamic surgical environments.

For further details about the Head tracking Mouse, please [visit the dedicated page](#).

Watch the Head Tracking Mouse video for more insights.