



## Detect intraoral camera kit with 19” monitor and white Easy Tray cart

Dental imaging kit including **DETECT wireless intraoral camera** with wide-angle lens and blue diagnostic light for detecting caries, tartar and plaque, **19” Visiomed Full HD medical monitor with protective glass** and **Easy Tray white stainless steel cart with monitor support**.

**Available on request**

**SKU:** OFES1004

**Categories:** Equipment Carts, Complete Camera Sets

### Description

---

#### Intraoral camera with monitor and Easy Tray cart

---

This **dental intraoral camera system** with wireless connection, monitor and Easy Tray cart is designed to improve patient communication and clinical documentation in dental practices.

The **DETECT WiFi intraoral camera** features a **high-definition wide-angle lens** and a dedicated **blue diagnostic light**, designed to facilitate the visual detection of **caries, tartar and plaque**. This technology allows clinicians to easily identify treatment areas and clearly show patients the condition of their teeth.

The system includes a **19” Visiomed Full HD medical monitor**, with **16:10 aspect ratio** and **protective glass**, ensuring sharp images and comfortable viewing during intraoral examinations.

All components are installed on an **electrified white stainless steel Easy Tray cart**, equipped with **monitor support** and designed to provide stability, hygiene and easy mobility between treatment rooms.

The result is a **mobile intraoral imaging system**, ready to use and ideal for improving visual diagnosis, patient communication and treatment presentation.

**Complete dental imaging system ready to use, without the need for additional equipment.**

#### Clinical applications

---

The wireless intraoral camera system is designed to support patient communication and clinical documentation in modern dental practices.

Displaying images on the monitor allows clinicians to show **caries, tartar, plaque and other dental conditions in real time**, helping patients better understand diagnoses and treatment plans.



The Easy Tray mobile cart allows the system to be easily moved between operatories, making it a practical solution for:

- dental practices
- dental clinics
- dental hygiene centres
- visual diagnosis and patient consultation

### **Package includes**

---

- **DETECT WiFi intraoral camera** with wide-angle lens and blue light for detecting caries, tartar and plaque – Tecnomed Italia code **AP295HD**
- **19" Visiomed Full HD monitor**, 16:10 format with protective glass – Tecnomed Italia code **MOT19SE**
- **Easy Tray electrified white stainless steel cart with monitor support** – Tecnomed Italia codes **OF1003BI.00** and **SVET9010L**