



Detect intraoral camera kit with 19” monitor and white Easy Tray cart

Dental imaging kit including **DETECT wireless intraoral camera** with wide-angle lens and blue diagnostic light for detecting caries, tartar and plaque, **19” Visiomed Full HD medical monitor with protective glass** and **Easy Tray white stainless steel cart with monitor support**.

Available

SKU: OFES1004

Categories: Complete Camera Sets, Equipment Carts

Description

Intraoral camera with monitor and Easy Tray cart

This **dental intraoral camera system** with wireless connection, monitor and Easy Tray cart is designed to improve patient communication and clinical documentation in dental practices.

The **DETECT WiFi intraoral camera** features a **high-definition wide-angle lens** and a dedicated **blue diagnostic light**, designed to facilitate the visual detection of **caries, tartar and plaque**. This technology allows clinicians to easily identify treatment areas and clearly show patients the condition of their teeth.

The system includes a **19” Visiomed Full HD medical monitor**, with **16:10 aspect ratio** and **protective glass**, ensuring sharp images and comfortable viewing during intraoral examinations.

All components are installed on an **electrified white stainless steel Easy Tray cart**, equipped with **monitor support** and designed to provide stability, hygiene and easy mobility between treatment rooms.

The result is a **mobile intraoral imaging system**, ready to use and ideal for improving visual diagnosis, patient communication and treatment presentation.

Complete dental imaging system ready to use, without the need for additional equipment.

Clinical applications

The wireless intraoral camera system is designed to support patient communication and clinical documentation in modern dental practices.

Displaying images on the monitor allows clinicians to show **caries, tartar, plaque and other dental conditions in real time**, helping patients better understand diagnoses and treatment plans.



The Easy Tray mobile cart allows the system to be easily moved between operatories, making it a practical solution for:

- dental practices
- dental clinics
- dental hygiene centres
- visual diagnosis and patient consultation

Package includes

- **DETECT WiFi intraoral camera** with wide-angle lens and blue light for detecting caries, tartar and plaque – Tecnomed Italia code **AP295HD**
- **19" Visiomed Full HD monitor**, 16:10 format with protective glass – Tecnomed Italia code **MOT19SE**
- **Easy Tray electrified white stainless steel cart with monitor support** – Tecnomed Italia codes **OF1003BI.00** and **SVET9010L**