



Ivory angle head, 1:1, prophylaxis use, compatible KaVo 31LR-L31

Ivory angle head, 1:1, prophylaxis use, Snap-On compatible KaVo 31LR-L31

Available

SKU: SL621P

Categories: KaVo - Ivory - Mk-dent

Description

The **Ivory 1:1 contra-angle head compatible with KaVo 31LR-L31** is designed to provide precision, smooth operation and comfort during oral hygiene and maintenance procedures.

When used with **Silk Line** contra-angles, it allows easier access to posterior quadrants, helping simplify daily clinical work and improve instrument handling during use.

The **rotostop** locking system allows quick and secure head replacement, making it easier to adapt to different transmission ratios and operating requirements.

Compatible with all **Ivory** and **KaVo®** contra-angles, it ensures excellent interchangeability and more efficient instrument management.

The **wear-resistant front bearing with integrated reinforcement bushing** provides stability, reduced vibration and long service life. The structure is designed to absorb stress during use, helping preserve rotor concentricity and maintain consistent performance.

Technical features Ivory 1:1 contra-angle head compatible with KaVo 31LR-L31

Indications	Standard
Ratio	2:1
Locking mechanism	Button
Head diameter (mm)	9,5
Head height (mm)	11,6
Bur diameter (mm)	2,334/2,350



Bur length (mm)	21.95
Type of lubrication	Oil
Autoclave resistant	Yes, 135°C
Packing size lxdxh (mm)	105 x 75 x 25 1
Packing weight appr. (g)	120
Certification	CE