




Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/ivory-angle-head-1-1-prophylaxis-use-compatible-kavo-31lr-l31/>



Ivory angle head, 1:1, prophylaxis use, compatible KaVo 31LR-L31

Ivory angle head, 1:1, prophylaxis use, Snap-On compatible KaVo 31LR-L31

218,14 €

Available

SKU: SL621P

Categories: Professional Equipment

Description

The Ivory angle head 1:1 prophylaxis use, compatible KaVo 31LR-L31.

The small size of the Ivory heads fitted to Silk Line prophylactic and endodontic contra-angles ensures easy access to preparation sites at the rear of all quadrants. Ivory heads are compatible with all Ivory and Kavo® contra-angles.

Thanks to the Rotostop locking system, heads can be removed and changed quickly and easily on all prophylactic line instruments to give the desired reduction ratio.

This high level of versatility permits quick and effective maintenance and lets you replace a damaged head with a new one without having to interrupt an operation if the need occurs. The front bearing features an integrated reinforcement bushing to absorb stress during strenuous operations and guarantee the concentricity of the rotor at all times.

Technical features Ivory angle head 1:1 prophylaxis use.

Indications	Standard
Ratio	2:1
Locking mechanism	Button
Head diameter (mm)	9,5
Head height (mm)	11,6
Bur diameter (mm)	2,334/2,350
Bur length (mm)	21.95
Type of lubrication	Oil
Autoclave resistant	Yes, 135°C
Packing size lxdxh (mm)	105 x 75 x 25 1
Packing weight appr. (g)	120
Certification	CE