



Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/ivory-angle-head-11-standard-use-compatible-kavo-80ld-l80/>



Ivory angle head, 1:1, standard use, compatible KaVo 80LD-L80

Ivory angle head, 1:1, standard use, FG shank 1.6 mm, compatible KaVo 80LD-L80

218,14 €

Available on request

SKU: SL680C

Categories: Professional Equipment

Description

Ivory angle head 1:1 compatible KaVo 80LD-L80.

The small size of the Ivory heads fitted to Silk Line prophylactic and endodontic contra-angles ensures easy access to preparation sites at the rear of all quadrants. Ivory heads are compatible with all Ivory and Kavo® contra-angles.

Thanks to the Rotostop locking system, heads can be removed and changed quickly and easily on all prophylactic line instruments to give the desired reduction ratio.

This high level of versatility permits quick and effective maintenance and lets you replace a damaged head with a new one without having to interrupt an operation if the need occurs. The front bearing features an integrated reinforcement bushing to absorb stress during strenuous operations and guarantee the concentricity of the rotor at all times.

Technical features Ivory angle head 1:1 compatible KaVo 80LD-L80

Indications	Standard
Ratio	1:1
Locking mechanism	Button
Head diameter (mm)	9,4
Head height (mm)	9,5
Bur diameter (mm)	2,334/2,350
Bur length (mm)	21.95
Type of lubrication	Oil
Autoclave resistant	Yes, 135°C
Max speed (RPM)	40.000
Packing weight appr. (g)	120

Certification	CE
---------------	----