



## Ivory 1:1 contra-angle head, KaVo compatible, NEW LINE CHANGE

Ivory angle head, 1:1, prophylaxis use, Snap-On compatible KaVo 31LR-L31. **BRAND NEW - PRODUCT LINE CHANGE ONLY 10 PIECES AVAILABLES.**

**Available**

**SKU:** OFOUT105

**Categories:** Outlet

### Description

Ivory 1:1 contra-angle head compatible KaVo 31LR-L31 – **BRAND NEW – PRODUCT LINE CHANGE ONLY 10 PIECES AVAILABLES.**

The head delivers precision, smooth performance and comfort during professional hygiene and maintenance procedures. When used with Silk Line contra-angles, it ensures easy access to posterior quadrants and preparation sites.

With the **rotostop locking system**, the head can be easily removed or replaced, adapting to different reduction ratios and clinical applications. It is **fully compatible with all Ivory and KaVo® contra-angles**, guaranteeing perfect interchangeability and seamless operation.

The **anti-wear front bearing** with integrated reinforcement bushing absorbs stress during use, maintaining rotor concentricity and extending the product's service life. Its robust structure minimises vibration and ensures long-term reliability.

#### Technical features Ivory 1:1 contra-angle head compatible KaVo 31LR-L31

Indications	Standard
Ratio	2:1
Locking mechanism	Button
Head diameter (mm)	9,5
Head height (mm)	11,6
Bur diameter (mm)	2,334/2,350
Bur length (mm)	21.95



Type of lubrication	Oil
Autoclave resistant	Yes, 135°C
Packing size lxdxh (mm)	105 x 75 x 25 1
Packing weight appr. (g)	120
Certification	CE