



## Ivory angle head, 3:1, surgical use, KaVo CL3 compatible

Ivory angle head, 3:1, surgery use, KaVo CL3 compatible

**Available on request**

**SKU:** SL631C

**Categories:** KaVo - Ivory - Mk-dent

### Description

The **Ivory 3:1 surgical contra-angle head compatible with KaVo CL3** is engineered for power, stability and precision during surgical procedures. Designed for **Ivory Silk Line contra-angles**, it ensures outstanding performance and simplified access to posterior quadrants thanks to its **compact dimensions**.

Fully compatible with **Ivory and KaVo® contra-angles**, it provides exceptional versatility. The **Rotostop locking system** allows quick and secure head replacement, adapting the reduction ratio as needed for each clinical application.

With its **3:1 reduction ratio**, this head delivers excellent torque control and smooth power transmission, making it ideal for surgical work and procedures requiring delicate handling. The **front anti-wear bearing** with **integrated reinforced bushing** absorbs operational stress, maintaining perfect rotor concentricity for quiet and long-lasting operation.

Its **modular construction** allows quick maintenance and easy replacement without interrupting workflow, ensuring efficiency and reliability in daily clinical practice.

Technical features **Ivory 3:1 surgical contra-angle head compatible with KaVo CL3:**

|                    |             |
|--------------------|-------------|
| Indications        | Surgery     |
| Ratio              | 3:1         |
| Locking mechanism  | Button      |
| Head diameter (mm) | 9,8         |
| Head height (mm)   | 10          |
| Bur diameter (mm)  | 2,334/2,350 |



|                          |            |
|--------------------------|------------|
| Bur length (mm)          | 21.95      |
| Type of lubrication      | Oil        |
| Autoclave resistant      | Yes, 135°C |
| Max speed (RPM)          | 10.000     |
| Packing weight appr. (g) | 120        |
| Certification            | CE         |