



Ivory angle head, 3:1, surgical use, KaVo CL3 compatible

Ivory contra-angle head, 3:1, surgery use, KaVo CL3 compatible

Estimated shipping in 15 days

SKU: SL631C

Categories: KaVo - Ivory - Mk-dent

Description

The **Ivory 3:1 surgical contra-angle head compatible with KaVo CL3** is designed to provide power, stability and precision during surgical procedures.

Ideal for **Ivory Silk Line** contra-angles, it delivers high performance and easier access to posterior quadrants thanks to its compact dimensions.

Compatible with all **Ivory** and **KaVo®** contra-angles, it offers excellent versatility. The **Rotostop** locking system allows the head to be removed and replaced quickly, adapting the reduction ratio according to the type of procedure and improving clinical productivity.

Thanks to the **3:1 transmission**, it provides optimal torque control and uniform speed reduction, making it ideal for surgical procedures and applications requiring precision and tactile sensitivity.

The **wear-resistant front bearing with integrated reinforcement bushing** reduces vibrations and ensures rotor concentricity, providing quiet and long-lasting operation.

The modular structure allows fast and targeted maintenance: when service is required, the head can be replaced without interrupting work, ensuring efficiency and operating continuity.

Technical features **Ivory 3:1 surgical contra-angle head compatible with KaVo CL3:**

Indications	Surgery
Ratio	3:1
Locking mechanism	Button
Head diameter (mm)	9,8



Head height (mm)	10
Bur diameter (mm)	2,334/2,350
Bur length (mm)	21.95
Type of lubrication	Oil
Autoclave resistant	Yes, 135°C
Max speed (RPM)	10.000
Packing weight appr. (g)	120
Certification	CE