



## Ivory Neo Line 1:1 blue contra-angle, spray head comp. KaVo

Ivory Neo Line 1:1 blue contra-angle for prophylaxis, Intramatic coupling, internal spray, Kavo compatible head

**Available**

**SKU:** SL411E

**Categories:** Contra-Angles handpieces

### Description

The **blue contra-angle Ivory Neo Line 1:1 blind compatible with Kavo, spray and removable head** is our direct-drive 1:1 contra-angle for prophylaxis, designed without optical fibre to deliver smooth, reliable and precise performance.

Like all instruments in the Ivory **Prophylaxis Line**, it combines high-quality engineering with advanced microtechnology to ensure quiet operation, exceptional control and outstanding balance in daily use.

Its **removable head** makes cleaning and maintenance easier while allowing the use of different head types — both standard and reduction — on the same handpiece, optimising investment and workflow. The compact head design improves access to sensitive or hard-to-reach areas for greater clinical accuracy.

The **high-precision transmission system** minimises mechanical tolerances, eliminating vibration, noise and power loss from the micromotor. This ensures a smooth and stable rotation under all operating conditions.

The Ivory Neo Line 1:1 blue contra-angle can be **thermo-disinfected and sterilised in an autoclave up to 135°C**, providing hygiene, reliability and durability.

Technical features Ivory Neo Line 1:1 blue contra-angle KaVo compatibility

Removable head	YES
Optical fiber	NO
Transmission	1: 1
CE	0476
Indications	Prophylaxis



Stop mechanism	Button
Type of Spray	Internal
Compatibility	Intramatic Universal
Rotation speed	40.000
Head Rotation	NO
Diameter of the Head	9.5mm
Head Height	11.6mm
Diameter of the Bur	2.334/2.350 mm
Length of the Bur	9-34mm
Weight	79gr
Autoclavable	YES max 135°C
Type of Lubrication	Oil