



## Ivory Neo Line 1:1 blue contra-angle, spray head comp. KaVo

Ivory Neo Line 1:1 blue contra-angle for prophylaxis, Intramatic coupling, internal spray, Kavo compatible head

**Shipping expected in approximately 15 days**

**SKU:** SL411E

**Categories:** Contra-Angles handpieces

### Description

The **blue Ivory Neo Line 1:1 blind blue contra-angle, KaVo compatible, with spray and removable head** is a direct-drive 1:1 contra-angle for prophylaxis, without fibre optics, designed to provide reliable performance, operating precision and excellent handling.

Manufactured with careful construction and advanced microtechnology, it offers smooth movement, precision and quiet operation during each procedure. The high-precision mechanical transmission maintains an extremely low tolerance threshold, minimising vibrations and energy loss.

The **removable head** makes cleaning and maintenance easier, helping to optimise costs thanks to the possibility of using interchangeable heads, both standard and reduction heads, with a single contra-angle body. The reduced size of the head also allows easier access to sensitive and hard-to-reach areas.

Operation is quiet, smooth and free from unwanted oscillations, providing stable and comfortable use during prophylaxis procedures.

The **blue Ivory Neo Line 1:1 blind contra-angle, KaVo compatible, with spray and removable head** can be thermally decontaminated and sterilised in an autoclave up to **135°C**, ensuring hygiene and long service life over time.

#### Technical features Ivory Neo Line 1:1 blue contra-angle KaVo compatibility

Removable head	YES
Optical fiber	NO
Transmission	1: 1
CE	0476



Indications	Prophylaxis
Stop mechanism	Button
Type of Spray	Internal
Compatibility	Intramatic Universal
Rotation speed	40.000
Head Rotation	NO
Diameter of the Head	9.5mm
Head Height	11.6mm
Diameter of the Bur	2.334/2.350 mm
Length of the Bur	9-34mm
Weight	79gr
Autoclavable	YES max 135°C
Type of Lubrication	Oil