




Printed on 28/05/2025

<https://shop.tecnomeditalia.com/en/product/ivory-silk-line-15-red-contra-angle-with-spray/>



## Ivory Silk Line 1:5 red contra-angle with spray

Ivory Silk Line 1:5 red multiplier contra-angle, Intramatic coupling, internal spray

**677,36 €**

**Available on request**

**SKU:** SL415C

**Categories:** Professional Equipment

### Description

The Ivory Silk Line 1:5 red contra-angle with spray is our 1: 5 multiplier contra-angle without optical fiber.

Perfect for all restorative operations and prosthetic implantation, it has been designed with a smaller head to offer you total freedom of movement.

Its Monoblock design facilitates a better cleaning, an even quieter operation and a high degree of manageability.

In addition, the anti-wear ceramic ball bearings guarantee a longer life cycle and a lower than average acoustic emission and its solid structure of the external parts in noble metal, guarantees for a long time its resistance to autoclave sterilization ... it is a piece that cannot be missing in your clinic!

Technical features Ivory Silk Line 1:5 red contra-angle with spray:

Removable head	NO
Optical fiber	NO
Transmission	1:5
Rotor bearings	In Ceramics
Indications	Prosthetic Restoration and Restoration
Stop mechanism	Button
Type of Spray	Internal
Compatibility	Intramatic Universal
Rotation speed	200.000
Head Rotation	NO
Diameter of the Head	9.9mm
Head Height	14.5mm
Diameter of the Bur	1.59-1.60mm

Length of the Bur	10-25mm
Weight	82gr
Autoclavable	YES
Type of Lubrication	Oil