



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/ivory-silk-oil-lubricating-oil-500ml-universal/>



Ivory Silk Oil lubricating oil, 500ml, universal

Ivory Silk Oil lubricating oil for handpieces, 500 ml spray, universal

17,16 €

Available

SKU: OL101S

Categories: Lubrification

Description

Ivory lubricating oil 500ml.

The Ivory range of lubricants is universally recognised for its unrivalled degreasing and cleaning power.

- elimination of overheating in moving parts;
- the prevention of deposits;
- lower consumption and running costs;
- extended instrument life;
- zero thermal stress and oxidation;
- better performance and fewer technical assistance callouts

just some of the benefits of using this oil.