

## Jaw, complete, diameter 45 mm, dental white



Fixing clamp for monitor support arm, designed for installation on dental unit pole with 45mm diameter.

**Available**

**SKU:** GA459010

**Categories:** Pole Clamps & Mounting Blocks

### Description

---

**GA459010 – Fixing clamp for monitor arm on 45mm diameter medical pole.**

#### Main features

- Vertical pole mounting
- Internal cable routing
- Clamping system with 2 × M8 screws

#### Technical features

- Compatible with Ø45mm poles
- Monitor arm insertion diameter: Ø35mm
- Dimensions (L × H): 111 × 90mm
- Material: aluminium
- Colour: dental white



- Certification: 2001/95/EC
- Warranty: 36 months

### **Applications**

- Dental and podiatry treatment units with Ø45mm lamp poles
- Equipped poles with Ø45mm diameter

### **Configuration note**

The fixing clamp for 45mm medical pole can be combined with Tecnomed monitor arms or equivalent monitor arms featuring a Ø35mm spigot.