



Jaw, complete, diameter 45 mm, dental white

Fixing block for monitor support arm, designed for installation on dental unit pole with 45mm diameter.

Available

SKU: BL459010

Categories: Pole Clamps & Mounting Blocks

Description

BL459010. fixing block for monitor arm on 45mm diameter pole.

Main features

- Vertical pole mounting
- Internal cable routing
- Clamping system with 2 × M8 screws

Technical features

- Compatible with Ø45mm poles
- Monitor arm insertion diameter: Ø35mm
- Dimensions (L × H): 52 × 90mm
- Material: aluminium
- Colour: dental white



- Certification: 2001/95/EC
- Warranty: 36 months

Applications

- Dental and podiatry treatment units with Ø45mm lamp poles
- Equipped poles with Ø45mm diameter

Configuration note

The fixing block can be combined with Tecnomed monitor arms or equivalent monitor arms featuring a Ø35mm spigot.