

Jaw, complete, diameter 50 mm, dental white



Fixing clamp for monitor support arm, designed for installation on dental unit pole with 50mm diameter.

Available from 10/09/2026

SKU: GA509010

Categories: Pole Clamps & Mounting Blocks

Description

GA509010 – fixing clamp for monitor arm on 50mm diameter medical pole.

Main features

- Vertical pole mounting
- Internal cable routing
- Clamping system with 2 × M8 screws

Technical features

- Compatible with Ø50mm poles
- Monitor arm insertion diameter: Ø35mm
- Dimensions (L × H): 116 × 90mm
- Material: aluminium
- Colour: dental white



- Certification: 2001/95/EC
- Warranty: 36 months

Applications

- Dental and podiatry treatment units with Ø50mm lamp poles
- Equipped poles with Ø50mm diameter

Configuration note

The fixing clamp for Ø 50mm medical pole can be combined with Tecnomed monitor arms or equivalent monitor arms featuring a Ø35mm spigot.