



## Jaw, complete, diameter 50 mm, dental white

Fixing block for monitor support arm, designed for installation on dental unit pole with 50 mm diameter.

**Available**

**SKU:** BL509010

**Categories:** Pole Clamps & Mounting Blocks

### Description

---

**BL509010 – fixing block for monitor arm on 50mm diameter pole.**

#### Main features

- Vertical pole mounting
- Internal cable routing
- Clamping system with 2 × M8 screws

#### Technical features

- Compatible with Ø50mm poles
- Monitor arm insertion diameter: Ø35mm
- Dimensions (L × H): 52 × 90mm
- Material: aluminium
- Colour: dental white



- Certification: 2001/95/EC
- Warranty: 36 months

### **Applications**

- Dental and podiatry treatment units with Ø50mm lamp poles
- Equipped poles with Ø50mm diameter

### **Configuration note**

The fixing block on 50mm pole can be combined with Tecnomed monitor arms or equivalent monitor arms featuring a Ø35mm spigot.