

## Jaw, complete, diameter 55 mm, dental white



Fixing clamp for monitor support arm, designed for installation on dental unit pole with 55mm diameter.

**Available**

**SKU:** GA559010

**Categories:** Pole Clamps & Mounting Blocks

### Description

---

**GA559010 – Fixing clamp for monitor arm on 55mm diameter poles**

#### Main features

- Vertical pole mounting
- Internal cable routing
- Clamping system with 2 × M8 screws

#### Technical features

- Compatible with Ø55mm poles
- Monitor arm insertion diameter: Ø35mm
- Dimensions (L × H): 121 × 90mm
- Material: aluminium
- Colour: dental white



- Certification: 2001/95/EC
- Warranty: 36 months

### **Applications**

- Dental and podiatry treatment units with Ø55mm lamp poles
- Equipped poles with Ø55mm diameter

### **Configuration note**

The fixing clamp can be combined with Tecnomed monitor arms or equivalent monitor arms featuring a Ø35mm connection pin.