



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/kit-intraoral-camera-ap295hd-equipment-tray-es1003bs/>



Kit intraoral camera AP295HD + equipment tray ES1003BS

Kit consisting of: Tecnomed Intraoral camera DETECT, equipment trolley ES1003BS, LCD monitor 19-inch, 1440*900 resolution, 16:10

1.851,81 €

Available on request

SKU: OFES1003

Categories: Multifunctional cart

Description

Kit intraoral camera DETECT + equipment tray SOFT

Kit consisting of:

- Tecnomed Intraoral camera DETECT;
- equipment trolley ES1003BS
- LCD monitor 19-inch, 1440*900 resolution, 16:10

– Tecnomed Intraoral camera DETECT

The DETECT intraoral camera has been designed to be an excellent diagnostic tool, offering detailed and high-quality images of the oral cavity. This device stands out for its ability to precisely identify occlusal and surface caries, plaque, and tartar. The ease of wireless image transfer makes the DETECT intraoral camera particularly practical and efficient.

Main features:

- Wide-angle lens with automatic focus.
- Specific blue light for detecting caries, tartar, and plaque.
- Ability to easily view and switch between one and four images on the screen.

Technical features:

- Functioning: autonomous and/or PC
- Sensor: 1/3 HD CMOS
- Resolution (pixel): 1920×1080

- Definition (Mp): 2.0
- Sensitivity/luminosity (lum): 6
- Automatic power off at rest: si
- Outputs: VGA/HDMI/WI-FI
- Handpiece Weight (gr): 45
- Handpiece connection: cable
- Cable lenght (mt): 2
- Illumination: 8 LED + LED DETECT BLUE
- Lens: fixed
- Perspective angle: 105°
- Autofocus (mm): 3-50
- Docking station: USB 2GB (up to 1000 images)
- Software: windows 10(32-64bit)- Windows 11
- Power supply: 12 VDC / 240 VAC
- Power supply for dental unit (DE5.0612): 12VDC/24 VAC
- Warranty: 24
- Certification: CE MDR 2017/745/EU

– LCD monitor 19-inch

Tecnomed Italia slim-line LCD monitor 19-inch feature high luminosity LCD screen. Tecnomed Italia monitor is fitted with a special outer glass to protect the LCD screen and ensure impermeability to cleaning solutions.

This monitor is designed for superb performance with both colour and greyscale images.

Technical features LCD monitor 19-inch

Flat Screen: 19" TFT LCD

Outer Glass: Glass

- Aspect Ratio: 16:9;
- Resolution: 1440 x 900;
- Power Supply: 12 VDC / 220 VAC;
- Mounting: VESA 100;
- Input: HDMI / DVI;
- Dimensions (W x L x H cm): 45 x 29 x 6
- Color: RAL9003;
- Warranty (months): 24;
- CE: Yes.

– Equipment trolley white, 3 shelves

The ES1003BS Dental Cart is a fundamental element for any medical or dental practice, laboratory or beauty center, thanks to its versatility and user-friendliness.

Main features Dental Cart SOFT 3 layers

- Made of metal, highly resistant;
- easy to clean and sanitize;
- Smooth movement;
- 3 layers (up 53x29,5 cm, middle 27,5x47,5 cm, down 54x34 cm);
- middle layer adjustable in height;

- dimensions 60x41x89cm.