



## Kit monitor + intraoral camera, polemount 45mm

Kit consisting of: LCD monitor 19-inch, 1440\*900 resolution, 16:10, Tecnomed Intraoral camera DETECT, 5 MP, blue light, HDMI & VGA output, WiFi, Arm for monitor, 42/45 mm pole mount, dental white.

**Available**

**SKU:** OF522HD

**Categories:** Complete Camera Sets

### Description

---

Kit monitor + intraoral camera DETECT + pole mount 45mm.

consisting of:

- LCD monitor 19-inch, 1440\*900 resolution, 16:10
- Tecnomed Intraoral camera DETECT, 5 MP, blue light, HDMI & VGA output, WiFi
- Arm for monitor, 42/45 mm pole mount, dental white

### Tecnomed Intraoral camera DETECT

---

The DETECT intraoral camera has been designed to be an excellent diagnostic tool, offering detailed and high-quality images of the oral cavity. This device stands out for its ability to precisely identify occlusal and surface caries, plaque, and tartar. The ease of wireless image transfer makes the DETECT intraoral camera particularly practical and efficient.

#### Main features:

- Wide-angle lens with automatic focus.
- Specific blue light for detecting caries, tartar, and plaque.
- Ability to easily view and switch between one and four images on the screen.



### **Technical features:**

- Functioning: autonomous and/or PC
- Sensor: 1/3 HD CMOS
- Resolution (pixel): 1920x1080
- Definition (Mp): 2.0
- Sensitivity/luminosity (lum): 6
- Automatic power off at rest: si
- Outputs: VGA/HDMI/WI-FI
- Handpiece Weight (gr): 45
- Handpiece connection: cable
- Cable length (mt): 2
- Illumination: 8 LED + LED DETECT BLUE
- Lens: fixed
- Perspective angle: 105°
- Autofocus (mm): 3-50
- Docking station: USB 2GB (up to 1000 images)
- Software: windows 10(32-64bit)- Windows 11
- Power supply: 12 VDC / 240 VAC
- Power supply for dental unit (DE5.0612): 12VDC/24 VAC
- Warranty: 24
- Certification: CE MDR 2017/745/EU

---

### **LCD monitor 19-inch**

Tecnomed Italia slim-line LCD monitor 19-inch feature high luminosity LCD screen. Tecnomed Italia monitor is fitted with a special outer glass to protect the LCD screen and ensure impermeability to cleaning solutions.

This monitor is designed for superb performance with both colour and greyscale images.

### **Technical features LCD monitor 19-inch:**

- Flat Screen: 19" TFT LCD;
- Outer Glass: Glass;
- Aspect Ratio: 16:9;
- Resolution: 1440 x 900;
- Power Supply: 12 VDC / 220 VAC;
- Mounting: VESA 100;
- Input: HDMI / DVI;
- Dimensions (W x L x H cm): 45 x 29 x 6;
- Color: RAL9003;
- Warranty (months): 24;
- CE: YES



## Arm for monitor, 42/45 mm pole mount, dental white

---