



Kit of 5 instrument trays, 415x300x20 mm

Kit including 5 pull-out instrument trays compatible with Easy Tray medical trolleys for workstation configuration.

Available on request

SKU: OF1003X.2

Categories: Stainless Steel Trays

Description

The **Easy Tray pull-out tray kit** 415x300x20mm includes 5 instrument trays for Easy Tray medical trolleys. It is designed to fully configure or reorganise the workstation, improving order, accessibility and workflow continuity.

Key features

- Kit including 5 pull-out trays
- Compatible with Easy Tray medical trolleys
- Suitable for clinical instruments and devices
- Easy-to-clean surfaces
- Single tray dimensions: 415x300x20mm