



Kit of 5 instrument trays, 465x315x20 mm

Kit including 5 pull-out instrument trays compatible with Easy Tray medical trolleys for workstation configuration.

Available

SKU: OF1003X.1

Categories: Stainless Steel Trays

Description

The **Easy Tray pull-out tray kit** 465x315x20mm includes 5 instrument trays designed for the organised placement of clinical instruments and devices on Easy Tray medical trolleys. It allows the configuration of a functional mobile workstation, ensuring practicality, order and efficient equipment management during clinical activities.

Key features

- Kit including 5 pull-out instrument trays
- Compatible with Easy Tray medical trolleys
- Support for clinical instruments and devices
- Easy-to-clean surfaces
- Single tray dimensions: 465x315x20mm