



## Kit Ultrasonic Cleaner 5 lt.

Composition Kit: 1 ultrasonic cleaner digiclean CLEA5L, 1 stand for 2 Pyrex glass CLA04, 1 narrow-mesh stainless steel basket CLA02N, 1 narrow-mesh cylindrical basket CLA02P, 2 universal Pyrex Glass CLA03.

**Available on request**

**SKU:** OFCLEA5L

**Categories:** Ultrasonic Cleaners

### Description

---

#### Kit Digiclean 5L Ultrasonic Cleaner – Fast, Deep and Reliable Cleaning in Stainless Steel

The **Digiclean 5L ultrasonic cleaner** delivers powerful and effective cleaning performance, combining high ultrasonic power with a disinfectant detergent. It removes dirt even from hard-to-reach areas, making it ideal for dental and medical instruments.

##### Durable Stainless Steel Construction

Made entirely of **scratch-resistant stainless steel**, the Digiclean 5L ensures long-term resistance to corrosion and wear. The **lid** and **drain tap** make daily operation and maintenance easy and efficient.

##### Intuitive Digital Controls

The **digital control panel** allows precise adjustment of **time (1–99 minutes)**, **temperature (0°C–80°C)** and **ultrasonic frequency**, ensuring a custom cleaning cycle for every application.

##### Complete Equipment Kit Digiclean 5L Ultrasonic Cleaner

Included in the package:

- Digiclean 5L ultrasonic cleaner
- Large and small fine-mesh stainless steel baskets



- Two 600ml Pyrex beakers
- Beaker support rack

Fully compliant with **EN ISO 13485:2016**, it offers maximum hygiene, durability and reliability for professional environments.