



Kit vasca ultrasuoni 5 litri

Composition Kit: 1 ultrasonic cleaner digiclean CLEA5L, 1 stand for 2 Pyrex glass CLA04, 1 narrow-mesh stainless steel basket CLA02N, 1 narrow-mesh cylindrical basket CLA02P, 2 universal Pyrex Glass CLA03.

Available

SKU: OFCLEA5L

Categories: Ultrasonic cleaners

Description

Kit ultrasonic cleaner 5lt.

The Ultrasonic washing is the most technologically advanced and cost-effective means of cleaning.

Ultrasonic cleaner Digiclean 5l generate an high ultrasonic power that associated with a disinfectant cleanser will ensure you a quick and effective cleaning action by removing dirt even in the hard- to-access parts.

Made of non-scratch stainless steel, it proves resistant over time against wear-and-tear corrosive agents.

An electronic generator operated by a microprocessor optimizes power performance while a display control panel allows the setting and the viewing of vibration times and liquid temperature values.

The 5 liters and 10 liters models are equipped with an exhaust washing liquid faucet so you can drain water quickly and easily.

Composition Kit ultrasonic cleaner 5lt.

- 1x ultrasonic cleaner digiclean CLEA5L;
- 1x stand for 2 Pyrex glass CLA04;
- 1x narrow-mesh stainless steel basket CLA02N;
- 1x narrow-mesh cylindrical basket CLA02P;
- 2x universal Pyrex Glass CLA03.