



## Laboratory turbine rotor for HCP11- HCP22S comp. Bien-Air TD783-TDS890

Dental laboratory turbine rotor compatible with Tecnomed HCP11 and HCP22S, Bien-Air TD783 and TDS890, supplied with ceramic torlon bearings.

**Available on request**

**SKU:** HBB837

**Categories:** TD783-TDS890

### Description

---

**Laboratory dental turbine rotor** compatible with Tecnomed HCP11 and HCP22S, Bien-Air TD783 and TDS890.

### Compatibility

- Tecnomed HCP11
- Tecnomed HCP22S
- Bien-Air TD783
- Bien-Air TDS890

### Technical features

- Ceramic torlon bearings
- Dynamically balanced rotor
- Serial number identification
- 12-month warranty
- Compliance with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1 and EN ISO 5349 standards