



## LCD monitor 19-inch, 1440\*900 resolution, 16:10

LCD monitor 19-inch, resolution 1440\*900, ratio 16:10 with special outer glass to protect the LCD screen, for dental and medical use.

**Available on request**

**SKU:** MOT19SE

**Categories:** Monitor-Pc

### Description

Tecnomed Italia slim-line LCD monitor 19-inch feature high luminosity LCD screen. Tecnomed Italia monitor is fitted with a special outer glass to protect the LCD screen and ensure impermeability to cleaning solutions.

This monitor is designed for superb performance with both colour and greyscale images.

Technical features LCD monitor 19-inch

Flat Screen	19" TFT LCD
Outer Glass	Glass
Aspect Ratio	16:9
Resolution	1440 x 900
Power Supply	12 VDC / 220 VAC
Mounting	VESA 100
Input	HDMI / DVI



Dimensions (W x L x H cm)	45 x 29 x 6
Color	RAL9003
Warranty (months)	24
CE	•