



## LCD monitor 19-inch, 1440\*900 resolution, 16:10 With medical power supply

LCD monitor 19-inch, resolution 1440\*900, ratio 16:10 with medical power supply and special outer glass to protect the LCD screen, for dental and medical use.

**Available on request**

**SKU:** MOT19MEN

**Categories:** Medical-Grade Monitors

### Description

#### Sanitisable Monitor VISIOMED – 19" High-Brightness LCD Display for Medical Environments

The **sanitisable monitor VISIOMED** delivers high image quality and hygienic safety in medical settings. Its 19" high-brightness LCD screen ensures precise and clear viewing during clinical procedures, even under challenging lighting conditions.

**Durable and easy to sanitise.**

A protective glass panel makes the screen resistant to scratches, water and impacts, allowing frequent disinfection with medical-grade solutions while maintaining visual clarity. Perfect for dental surgeries, clinics and laboratories.

**Superior visual performance.**

The 19" LED display in 16:10 format with Full HD resolution offers sharp, true-to-life images. VGA and HDMI inputs guarantee full compatibility with most medical and dental systems.

**Flexible and reliable installation.**

Equipped with a VESA 100 mount and 12VDC power supply, the VISIOMED monitor integrates seamlessly into wall mounts, arms or mobile carts, ensuring stability and adaptability in any professional setting.

#### Technical features LCD monitor 19-inch with medical power supply

Flat Screen	19" TFT LCD
-------------	-------------



Outer Glass	Glass
Aspect Ratio	16:9
Resolution	1440 x 900
Power Supply	12 VDC / 220 VAC
Medical mains adapter	IEC 60601
Mounting	VESA 100
Input	HDMI / DVI
Dimensions (W x L x H cm)	45 x 29 x 6
Color	RAL9003
Warranty (months)	24
CE	•