



## LCD monitor 19-inch, 1440\*900 resolution, 16:10 With medical power supply

LCD monitor 19-inch, resolution 1440\*900, ratio 16:10 with medical power supply and special outer glass to protect the LCD screen, for dental and medical use.

**Available on request**

**SKU:** MOT19MEN

**Categories:** Monitor-Pc

### Description

Tecnomed Italia slim-line LCD monitor 19-inch with medical power supply feature high luminosity LCD screen. Tecnomed Italia monitor is fitted with a special outer glass to protect the LCD screen and ensure impermeability to cleaning solutions.

This monitor is designed for superb performance with both colour and greyscale images.

Technical features LCD monitor 19-inch with medical power supply

Flat Screen	19" TFT LCD
Outer Glass	Glass
Aspect Ratio	16:9
Resolution	1440 x 900
Power Supply	12 VDC / 220 VAC
Medical mains adapter	IEC 60601
Mounting	VESA 100



Input	HDMI / DVI
Dimensions (W x L x H cm)	45 x 29 x 6
Color	RAL9003
Warranty (months)	24
CE	•