



LCD monitor 19-inch, 1440*900 resolution,16:10 With medical power supply

LCD monitor 19-inch, resolution 1440*900, ratio 16:10 with medical power supply and special outer glass to protect the LCD screen, for dental and medical use.

Available on request

SKU: MOT19MEN

Categories: Medical-Grade Monitors

Description

Sanitisable Monitor VISIOMED – 19” High-Brightness LCD Display for Medical Environments

The **sanitisable monitor VISIOMED** delivers high image quality and hygienic safety in medical settings. Its 19” high-brightness LCD screen ensures precise and clear viewing during clinical procedures, even under challenging lighting conditions.

Durable and easy to sanitise.

A protective glass panel makes the screen resistant to scratches, water and impacts, allowing frequent disinfection with medical-grade solutions while maintaining visual clarity. Perfect for dental surgeries, clinics and laboratories.

Superior visual performance.

The 19” LED display in 16:10 format with Full HD resolution offers sharp, true-to-life images. VGA and HDMI inputs guarantee full compatibility with most medical and dental systems.

Flexible and reliable installation.

Equipped with a VESA 100 mount and 12VDC power supply, the VISIOMED monitor integrates seamlessly into wall mounts, arms or mobile carts, ensuring stability and adaptability in any professional setting.

Technical features LCD monitor 19-inch with medical power supply

Flat Screen	19? TFT LCD
-------------	-------------



Outer Glass	Glass
Aspect Ratio	16:9
Resolution	1440 x 900
Power Supply	12 VDC / 220 VAC
Medical mains adapter	IEC 60601
Mounting	VESA 100
Input	HDMI / DVI
Dimensions (W x L x H cm)	45 x 29 x 6
Color	RAL9003
Warranty (months)	24
CE	•