



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/led-b-woodpecker-curing-light-1100-mw-cm2-cordless/>



LED-B Woodpecker curing light, 1100 mW/cm², cordless

LED-B Woodpecker curing light, 1100 mW/cm², cordless, 4 curing times

130,74 €

Available from 20/06/2025

SKU: PP630L

Categories: Cordless

Description

LED-B Woodpecker Curing Light: Quick and Precise Curing for High-Quality Dental Restorations.

The LED-B Woodpecker is a cordless curing light that ensures fast and even curing of dental composites, thanks to its advanced technology that maintains consistent light intensity without overheating. Ideal for dental practices looking for efficiency and reliability in their treatments.

Key Features:

- **Ergonomic and Comfortable Design:** The handpiece-like design provides excellent ergonomics, allowing for prolonged use without fatigue. The comfortable grip enhances precision during work.
- **Intuitive Control:** With an easy-to-use control panel, you can effortlessly select from four curing times: 5, 10, 15, and 20 seconds.
- **Safety and Practicality:** The black glass fiber optic tip is autoclavable at 135°C, ensuring maximum hygiene. The memory-free battery guarantees long-lasting performance and reliability.
- **Efficiency Check:** The included tester allows you to verify the correct functionality of the lamp for optimal curing results.

Technical Specifications:

- Light intensity: 1100 mW/cm²
- Curing times: 5, 10, 15, 20 seconds
- Tip autoclavable up to 135°C
- Power supply: memory-free battery
- Package includes: Handpiece, base, 8 mm tip, anti-glare shield, power adapter, tester

The LED-B Woodpecker is the ideal choice for professionals seeking high performance and reliability in their daily work.

Shop DTE-Woodpecker products at Tecnomed Italia and enhance your daily operations with the most advanced curing lights available, all supported by the reliable after-sales service from Tecnomed Italia.

Technical features LED-B Woodpecker curing light 1100 mW/cm2.

Light Source	5 W LED
Wavelength	420-490 nm
Light intensity	1000 mW/cm2 +/-10%
Handpiece dimensions	31 x 34 x 265 mm
Handpiece weight	110 gr.
Warranty	24 months