



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/lucy-wax-trolley-white-5-shelf-380x270x950-mm-lxpxh/>



Lucy wax trolley, white, 5 shelf, 380x270x950 mm (lxpxh)

Lucy wax trolley, 5 shelf, 380x270x950mm(lxpxh), compact size, metal structure on castors, 3 side shelves that can be closed.

407,90 €

Available on request

SKU: ES1004BC

Categories: Multifunctional cart

Description

Lucy wax trolley 5 shelf is sober and elegant design in a compact format. Practical and easy to handle. This is the ideal work tool.

The metal structure on castors has been designed to be a valid support in hair removal treatments.

It can in fact house a wax heater, rollers, and all of the other equipment needed for a professional job.

The Lucy wax trolley 5 shelf is comprised of a central sliding top, a fixed shelf on top, and three side shelves that can be closed; it also has a flip top container for a rubbish bag. The trolley can also be used for traditional treatments.

Dimensions: 38 x 27 cm

Height: 95 cm

Weight: 12 kg