



Lucy wax trolley, white, 5 shelf, 380x270x950 mm (lpxh)

Lucy wax trolley, 5 shelf, 380x270x950mm(lpxh), compact size, metal structure on castors, 3 side shelves that can be closed.

Estimated shipping in 21 days

SKU: ES1004BC

Categories: Beauty & Spa Treatment Furniture

Description

The **Lucy white wax trolley with 5 shelves** is a **practical and functional workstation** designed for **beauty salons and professional waxing treatments**.

Its **compact size** (380x270x950 mm) and **metal frame with wheels** make it easy to move and perfectly suited for small treatment rooms.

With a **clean and elegant design**, this trolley provides convenient organisation of waxing equipment and accessories.

It includes a **sliding central tray**, a **fixed top shelf**, and **three folding side trays**, allowing full flexibility during treatment preparation and application.

A **built-in tilting waste bin** ensures cleanliness and convenience, while the **five spacious shelves** can accommodate **wax heaters, applicators, tools and products** essential for professional use.

Constructed from **durable, easy-to-clean materials**, the **Lucy wax trolley 5 shelves** guarantees reliability and stability.

Weighing only **12 kg**, it combines **lightness and strength**, making it ideal not only for waxing procedures but also as a **multi-purpose beauty trolley** for aesthetic and wellness treatments.