



Professional waxing trolley with 5 shelves

Compact and practical professional waxing trolley with 5 shelves, metal frame on wheels, sliding central shelf, three foldable side shelves and concealed tilting waste bin holder. Made in Italy

Shipping expected in approximately 21 days

SKU: ES1004BC

Categories: Beauty & Spa Treatment Furniture

Description

This professional waxing trolley with 5 shelves is designed for beauty salons, spas and professional environments dedicated to waxing, hair removal and wellness treatments.

The metal frame on wheels makes it stable and easy to move during daily work, even in smaller treatment rooms. Its compact dimensions allow practical organisation of wax heaters, products and professional accessories.

The trolley features a fixed top shelf, a sliding central shelf and three foldable side shelves, useful for extending the working surface only when needed. The concealed tilting waste bin holder helps keep the workstation tidy and hygienic during treatments.

Thanks to its functional layout and compact structure, it can also be used as a multi-purpose trolley for beauty and treatment rooms. Made in Italy

TECHNICAL SPECIFICATIONS

- Top shelf: fixed
- Central shelf: sliding
- Side shelves: 3 foldable shelves
- Waste bin holder: concealed tilting system
- Colour: white
- Dimensions: 38 x 27 cm
- Height: 95 cm
- Weight: 12 kg
- Production: Made in Italy