



Massage table with adjustable backrest and face port, colour white

Massage table with face port, dimensions 1880x700x h790mm, sky colour white, adjustable backrest

Available on request

SKU: ES1007BL

Categories: Armchairs and massage tables

Description

Massage table with adjustable backrest.

Massage table with simple, stylish metal framework, with reclining headrest.

Ideal for beauty treatments and specialist exams.

Provided with paper roll holder and face-hole.

Optional: storage shelf.

Dimensions: 188 cm x 70 cm

Height: 79 cm

Weight: 53 Kg