

## Massage table with adjustable backrest and face port, colour white



Massage table with face port, dimensions 1880x700x h790mm, sky colour white, adjustable backrest

**Available on request**

**SKU:** ES1007BL

**Categories:** Beauty & Spa Treatment Furniture

### Description

---

The **white massage table with adjustable backrest and face hole** is a **versatile and functional solution** designed for **beauty salons, spas, medical clinics and physiotherapy studios**.

Its **simple yet sturdy metal frame** combines reliability, comfort and easy maintenance, ensuring long-lasting professional performance.

The **soft white sky upholstery** provides a hygienic, easy-to-clean surface that guarantees a pleasant and comfortable experience for the client. The **adjustable backrest** allows different reclining positions depending on the treatment, offering excellent support during massage, facial or medical procedures.

The **integrated face hole** ensures a relaxed and ergonomic posture when lying face down, making treatments more comfortable for both the operator and the client.

The table is supplied with a **paper roll holder**, ideal for quick preparation between sessions and maintaining high hygiene standards.

An **optional storage shelf** can be added to keep tools and products within reach, optimising organisation and workflow.

With its **dimensions of 1880x700x790 mm** and a **weight of 53 kg**, this massage table with adjustable backrest fits perfectly into professional environments that require **clean design, stability and functionality**.