



## Massage table with adjustable backrest, colour white

Massage table, dimensions 1820x620x h750mm, white sky, adjustable backrest

**Available on request**

**SKU:** ES1006BL

**Categories:** Beauty & Spa Treatment Furniture

### Description

---

The **white massage table with adjustable backrest** is a strong and practical choice for **beauty salons, spas and physiotherapy clinics**.

Available in **single- or double-section versions**, it offers maximum versatility for different treatments and working positions.

Its **powder-coated metal frame** ensures strength and stability, while the **washable eco-leather upholstery** provides comfort, hygiene and durability.

The **ergonomic backrest adjustment** allows operators to adapt the position according to each treatment, making it ideal for massage, beauty and therapeutic sessions.

In the **double-section model**, the lower leg area can be raised to improve **lymphatic drainage** or provide enhanced relaxation during treatment.

The table comes equipped with a **paper roll holder**, and an **optional storage shelf** can be added to keep accessories and tools within reach.

Compact yet stable, with dimensions of **182x62x75 cm** and a **weight of 42 kg**, this **massage table with adjustable backrest** fits perfectly into any professional environment dedicated to wellness, aesthetics and physiotherapy.