



Printed on 24/05/2025

<https://shop.tecnomeditalia.com/en/product/massage-table-with-adjustable-backrest-colour-white/>



Massage table with adjustable backrest, colour white

Massage table, dimensions 1820x620x h750mm, white sky, adjustable backrest

825,28 €

Available on request

SKU: ES1006BL

Categories: Armchairs and massage tables

Description

Massage table adjustable backrest.

Sturdy massage table with 1 or 2 jointed sections, ideal for beauty or physiotherapy treatment centres.

The whole framework is in painted metal and upholstered in washable faux leather.

It is provided with paper roll holder and can also be fitted with a practical storage compartment.

The twojointed section model offers the possibility to raise the lower leg section for lymph drainage or to add to the relaxation of different treatments.

Technical features Massage table adjustable backrest.

Optional: storage shelf.

Dimensions: 182 cm x 62 cm

Height: 75 cm

Weight: 42 Kg