

## Massage table with adjustable backrest, colour white



Massage table, dimensions 1820x620x h750mm, white sky, adjustable backrest

**Available on request**

**SKU:** ES1006BL

**Categories:** Armchairs and massage tables

### Description

---

Massage table adjustable backrest.

Sturdy massage table with 1 or 2 jointed sections, ideal for beauty or physiotherapy treatment centres.

The whole framework is in painted metal and upholstered in washable faux leather.

It is provided with paper roll holder and can also be fitted with a practical storage compartment.

The twojointed section model offers the possibility to raise the lower leg section for lymph drainage or to add to the relaxation of different treatments.

Technical features Massage table adjustable backrest.

**Optional:** storage shelf.

**Dimensions:** 182 cm x 62 cm

**Height:** 75 cm

**Weight:** 42 Kg