



Dental Mayo table in stainless steel with hydraulic system T-shaped base

Stainless steel AISI 304 **Mayo table** Made in Italy with **hydraulic foot pedal height adjustment**, T-base, **locking castors** and **removable tray**. Designed for stability, hygiene and ergonomic use in clinical environments.

Available on request

SKU: DE1006X

Categories: Medical & Dental Furniture

Description

Stainless steel Mayo table with hydraulic height adjustment

The **Mayo table in stainless steel AISI 304** is designed to provide a stable and adjustable support surface for instruments and materials during clinical procedures.

The **hydraulic foot pedal system** allows height adjustment up to **129 cm**, ensuring ergonomic use while keeping hands free during operations.

AISI 304 stainless steel structure and T-base

The structure in **AISI 304 stainless steel with scotch brite finish** ensures high corrosion resistance, durability and maximum hygiene in healthcare environments.

The **T-shaped base** provides stability and easy positioning even in limited spaces, while the **locking castors** allow smooth movement and secure positioning during use.

Sliding and removable tray

The table is equipped with a **removable tray**, allowing flexible instrument management during procedures.

This configuration improves accessibility and workflow efficiency, making the table suitable for professional use across different medical environments.



Technical features

- Structure in **AISI 304 stainless steel**
- **Scotch brite finish**
- **Hydraulic foot pedal height adjustment**
- Height adjustable up to **129 cm**
- **T-shaped base**
- **Locking castors**
- **Removable tray**
- Easy to clean and sanitise
- **Made in Italy**

Applications

- Dental clinics
- Medical clinics
- Veterinary clinics
- Operating rooms