



## Metasys H1 disinfectant for special areas, kit

Metasys H1 disinfectant for special areas, kit 4 x 500 ml pouches, for suction systems.

**Available on request**

**SKU:** 60010021

**Categories:** Special Applications

### Description

---

Metasys H1 disinfectant for special areas.

Decontamination concentrate for the METASYS automatic hygiene system H1, designed for the continuous cleaning of suction systems.

#### Key Features:

- **Continuous cleaning:** ensures the ongoing cleanliness of the entire suction system.
- **Active defoaming:** helps maintain system efficiency.
- **Certification:** CE 0044.

**Modern active agent combination:** quaternary ammonium compounds, special tensides, and enzymes.

**Packaging:** Kit with 4 pouches of 500 ml each.