



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/metasys-h1-disinfectant-for-special-areas-kit/>



Metasys H1 disinfectant for special areas, kit

Metasys H1 disinfectant for special areas, kit 4 x 500 ml pouches, for suction systems

157,84 €

Available on request

SKU: 60010021

Categories: Special areas

Description

Metasys H1 disinfectant for special areas.

Decontamination concentrate for the METASYS automatic hygiene system H1, designed for the continuous cleaning of suction systems.

Key Features:

- **Continuous cleaning:** ensures the ongoing cleanliness of the entire suction system.
- **Active defoaming:** helps maintain system efficiency.
- **Certification:** CE 0044.

Modern active agent combination: quaternary ammonium compounds, special tensides, and enzymes.

Packaging: Kit with 4 pouches of 500 ml each.