



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-ik-instrument-disinfectant-for-castellini-turbines/>



Metasys IK instrument disinfectant, for Castellini turbines

Metasys IK instrument disinfectant for Castellini turbines, 2 spray bottles 200 ml

91,37 €

Available on request

SKU: 60030187

Categories: Instruments

Description

Metasys IK instrument disinfectant for Castellini turbines.

Alcohol-based spray disinfectant, ideal for the rapid disinfection of water and air conduits in handpieces, contra-angles, and turbines. Its advanced formulation ensures effective disinfection in short times, providing enhanced safety through fixed adapters.

Key Features:

- **Rapid Efficiency:** Short action time, ideal for optimizing operational efficiency.
- **Optimized Safety:** Includes universal adapter for Castellini turbine.s
- **Extended Duration:** A single spray can provides approximately 300 applications (press for 2 seconds each time).
- **Active Ingredient Combination:** Advanced formula containing 1-propanol, 2-propanol, 2-butanol, and quaternary ammonium compounds.

Note: This product does not replace the lubrication recommended by the manufacturer.

Certification: CE 0044

Packaging Metasys IK instrument disinfectant for Castellini turbines: 2 x 200 ml spray bottles with adapters.