



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/metasys-ik-instrument-disinfectant-for-kavo-and-mk-dent-turbines/>



Metasys IK instrument disinfectant, for KaVo and MK-dent turbines

Metasys IK instrument disinfectant, for KaVo and MK-Dent turbines, 2 x 200 ml spray bottles with adapters

91,37 €

Available on request

SKU: 60030185

Categories: Instruments

Description

Metasys IK instrument disinfectant for KaVo and MK-Dent turbines, 2 x 200 ml spray bottles with adapters.

Alcohol-based spray disinfectant, ideal for the rapid disinfection of water and air conduits in handpieces, contra-angles, and turbines. Its advanced formulation ensures effective disinfection in short times, providing enhanced safety through fixed adapters.

Key Features:

- **Rapid Efficiency:** Short action time, ideal for optimizing operational efficiency.
- **Optimized Safety:** universal adapter compatible with KaVo and MK-Dent turbines.
- **Extended Duration:** A single spray can provides approximately 300 applications (press for 2 seconds each time).
- **Active Ingredient Combination:** Advanced formula containing 1-propanol, 2-propanol, 2-butanol, and quaternary ammonium compounds.

Note: This product does not replace the lubrication recommended by the manufacturer.

Certification: CE 0044

Packaging Metasys IK instrument disinfectant for KaVo and MK-Dent turbines: 2 x 200 ml spray bottles with adapters.